


Column : AppliChrom SugarSep-H V; 300mm x 8mm; P/N SASHV103008 Detection : RI
 Mobile Phase : H2O / 5mM H2SO4 Temperature : 40C
 Flow Rate : 0.4 Pressure : -
 Note : Sample Test Mix / Organic Acids



Clarity Lite - 16.05.2023 09_39_54

	Reten. Time [min]	Compound Name	Efficiency [th.p]	Eff/l [t.p./m]	Resolution [-]
1	14,52	Citric Acid	6454	21513	
2	16,22	Malic Acid	7957	26522	2,35
3	18,97	Succinic Acid	9475	31582	3,66
4	19,75	Lactic Acid	7448	24826	0,93
5	21,48	Formic Acid	10116	33721	1,96
6	23,14	Acetic Acid	10071	33570	1,87
7	25,07	Levulinic Acid	6291	20969	1,77
8	31,46	Iso-Butyric Acid	5601	18670	4,35
9	50,22	Valeric Acid	2587	8622	6,68